

GT 110


TESTA PER REFILARE E SMUSSARE

Cutterhead for trimming and chamfering
Messerkopf für beschneiden und kanten
Cabezal para ranurar y perfilar




CORPO IN LEGA LEGGERA, NON USARE SODA CAUSTICA

Light alloy body, do not use caustic soda
Körper aus Leichtmetall-Legierung, kein Ätznatron verwenden
Cuerpo en aleación ligera, no utilizar sosa cáustica



ESEMPI DI LAVORAZIONE

Types of processing
Verarbeitungsbeispiele
Ejemplos de procesamiento



CODICE CODE	MATERIALE MATERIAL	D mm	B mm	d mm	DKN mm	Z
GT11001	LEGA LEGG.	60	20	14	2 - 5x2,5	4
GT11002	LEGA LEGG.	70	20	16	2 - 5x2,5	6
GT11003	LEGA LEGG.	70	20	20	2 - 6x3	6
GT11004	LEGA LEGG.	80	20	16	2 - 5x2,5	6
GT11005	LEGA LEGG.	80	20	20	2 - 6x3	6

DISTINTA BASE/Bill of Materials

DESCRIZIONE DESCRIPTION	DIMENSIONI mm	CODICE CODE	01	02	03	04	05
Coltello/Knife	20x12x1,5	CGZ207	4	6	6	6	6
Cuneo/Wedge	18x12x5	CU027	4	6	6	6	6
Vite per Cuneo/Wedge's Screw	M8x10	VT015	4	6	6	6	6
Chiave per Cuneo/Wedge's Key	Esagonale 4	CH004	1	1	1	1	1